



GTIN: 7640178010413

Medi-Glaze

Dental Glazing resin

INTRODUCTION

Medi Glaze is a single-component, light-cured glazing agent enabling simultaneous treatment of both dentin and enamel using one liquid. It provides a clear and smooth surface.

II. INDICATIONS

Medi Glaze is indicated for the following restorations:
Surface treatment and glazing agent of porcelain, ceramics, hybrid ceramics or composite resin.
Direct and indirect composite resins.

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECTS

The oral mucosa membrane may turn slightly white when exposed to Medi Glaze, due to the coagulation of protein

V. INCOMPATIBILITIES

- 1- Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- 2- Do not use hemostatics containing ferric compounds, since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva, due to remaining ferric ions.
- 3- When using hemostatics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherent surface. Failure to do so might weaken the bond strength to the tooth structure.

VI. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue the use of the product, remove the product and seek medical attention.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye
4. If the product comes in contact with human body tissues, take the following actions: Immediately wash the eye with copious amounts of water and consult a physician.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the dental curing light when curing the product.
7. Avoid using the same product dispensed into a well of the dispensing dish and the same applicator brush for different patients to prevent cross infection. Discard the applicator brush after use.
8. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.

2. Handling and manipulation precautions

1. The product should not be used for purposes other than specified in [II.INDICATIONS].
2. The use of the product is restricted to licensed dentist.
3. The product contains ethanol, a flammable substance. Do not use it near an open flame.
4. Treat appropriately in a cavity close to the pulp or accidental pulp exposure with a pulp capping agent.
5. To prevent poor performance and handling characteristic, observe the specified light curing times and other handling requirements.
6. Clean the cavity sufficiently to prevent poor bonding. If the adherent surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
7. Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 7 minutes after dispensing. As the volatile ethanol which is contained in the **Medi Glaze** evaporates, the viscosity increases, thereby making it difficult to apply.
8. During the application of the **Medi Glaze** to the adherent surface and the duration of leaving it in place for 10 seconds after the application, move the lighting-spot out of the mouth or turn off the light to prevent the applied **Medi Glaze** from being exposed to the operating light because the **Medi Glaze** may gelate when left under such a light. In addition, after treatment for 10 seconds, dry sufficiently as soon as possible by blowing mild air.
9. **Medi Glaze** contains ethanol and water. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move; otherwise, the adhesion effect will be impaired. In order to dry sufficiently, adjust the air pressure according to the shape and size of the cavity and the prosthetic appliance. Use a vacuum aspirator to prevent the bond liquid from scattering.

10. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with the product.

11. If there is an adherent surface including uncut enamel, apply **Medi Etch GEL** to the uncut enamel, leave it in place for 10 seconds, wash with water, and then dry. After then, apply the product to an adherent surface.

12. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. The product must be refrigerated (2-8°C) when not in use, and should be brought to room temperature for more than 15 minutes before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature; otherwise, the excessive amount of liquid might be dispensed or the liquid might ooze after use.
3. The container should be tightly capped immediately after use to reduce the evaporation of the volatile solvent (ethanol contained in the product) . If the liquid does not flow easily from the nozzle, the nozzle may have plugged up. Do not dispense the liquid forcibly from the plugged container.
4. If the liquid may not flow easily, since it has not been used for a long term, shake the container before use.
5. Keep away from extreme heat, direct sunlight or a flame.
6. The product must be stored in proper places where only dental practitioners can access it.

VII. COMPONENTS

Principal ingredients:

- Bisphenol A diglycidylmethacrylate
- 2-Hydroxyethyl methacrylate
- Colloidal silica
- dl-Camphorquinone
- Ethanol

VIII. CLINICAL PROCEDURES

Application of **Medi Glaze**

1. Dispense the necessary amount of Medi Glaze into a well of the dispensing dish immediately before application.
2. Apply **Medi Glaze** to the entire cavity wall with the applicator brush. Leave it in place for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
3. Dry the entire cavity wall sufficiently by blowing mild air for more than 5 seconds until the bond does not move. Use a vacuum aspirator to prevent the bond liquid from scattering.
4. Light-cure **Medi Glaze** for 10 seconds with a dental curing unit (see table "Dental curing unit")

Table: Dental curing unit

Light power	500mW/cm ²	800mW/cm ²
Medi Glaze	40 sec.	20 sec.

Evaluated according to ISO 10650-1 & ISO 10650-2.

4. Polish with silicon rubber points or polishing discs in the usual manner.

[Caution]

[WARRANTY]

SwissMediTec GmbH will replace any product that is proven to be defective. **SwissMediTec GmbH** does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.



SwissMediTec GmbH
Businesspark
Bahnhofstrasse 42
6162 Entlebuch
Switzerland
Tel.: +41 41 481 00 86
+41 79 482 32 98
www.swiss-meditec.ch
Email: info@swiss-meditec.ch